

# SureTact<sup>®</sup> G3

## Sectional Matrix System

Predictability  
begins here





# SureTact® G3

## Sectional Matrix System

### Producing superior Class II Restorations

The SureTact® G3 sectional matrix system is a pioneering system comprising of Nickel-Titanium rings, anatomically contoured matrix bands, and adaptable self-guiding wedges that delivers dependable, snug contacts and restorations that precisely mimic the natural anatomy of the tooth.

*"Most clinicians find contact creation to be the most challenging part of a Class II restoration."*

The SureTact® G3 matrix system produces enhanced results, even in cases of extensive cavities and absent cusps. The shaped design of the matrices and wedges minimizes the need for restoration finishing. Sectional Matrix System has been trusted by tens of thousands of dental practitioners in clinics and educational institutions worldwide. The SureTact® G3 system sets the new standard for superior dentistry.



*Inspired by dentists to empower dentists*

# Predictable outcomes with SureTact® G3 Sectional Matrix System

Smart-Loop™ for easy insertion and removal of matrix


Strong ELAMAX Rings for outstanding spring strength and memory

Greater curvature: matrix wraps around tooth

Double U-shaped with ATR tines to accommodate the wedge


Gingival extension to prevent gaps in gingival-axial corner

Pronounced marginal ridge for natural-looking anatomy



## 4 simple steps to amazing Class II's

Refer to the Directions for Use for full instructions and reprocessing information.




1  
Place WedgeShield® prior to preparation to pre-separate and protect the adjacent tooth. Prepare tooth with conservative Class II preparation. Once the preparation is complete, remove the WedgeShield® and place your matrix.



2  
Place SureTact G3™ sectional matrix band. Select a matrix band that most closely approximates the occlusogingival height of the tooth. The Smart-Loop™ should be oriented toward the occlusal margin and may be bent over the adjacent tooth after placement.



3  
Insert SureTact G3™ wedge. Keep a finger on the matrix band to prevent the wedge from dislodging the band during insertion.



4  
Place the SureTact® G3 Ring on top of the wedge to provide separation

# Getting started

Each element within the SureTact® G3 Sectional Matrix System has been fine-tuned to simplify each stage of the process. Arrange a demonstration to witness for yourself how the components of the SureTact® G3 compare to your existing system.



## SureTact™ G3 Sectional Matrix System Intro Kit

100 Matrices – 25 each size: 3.5 mm, 4.5 mm, 5.5 mm, 6.5 mm  
75 Wedges – 25 each size: Small, Medium, Large  
2 Universal Ring; 1 Forceps; 1 Pin Tweezers

**M1001** Introductory Kit



## Universal Ring (Blue)

- ▶ Nickel-titanium for excellent holding force with ELAMAX technology
- ▶ ONE ring for all design to provide simple operation
- ▶ Stackable for multiple restorations
- ▶ Double U-shaped with ATR tines to provide best contact
- ▶ Advanced Thermal Resistant material was developed in Belgium

## Ring Refills

**M1002** Universal (2)



## MirrorCoat™ Matrices (4 sizes)

- ▶ Ultra-thin, MirrorCoat™ makes matrix removal easier and reduces light reflection, improving depth of field for better visibility
- ▶ Greater curvature for anatomical precision and less interference with ring placement
- ▶ Smart-Loop™ for precise insertion and removal
- ▶ XtraTite™ technology combines specially processed stainless steel with a proprietary matrix band forming process

## Matrix Refills

**M1004** Size 3.5 mm (100)

**M1005** Size 4.5 mm (100)

**M1006** Size 5.5 mm (100)

**M1007** Size 6.5 mm (100)



## Wedges (3 sizes)

- ▶ Wave-shaped wings compress on entry and flare upon exit for easy placement and a true fit and seal
- ▶ Hollow underside allows placement of second wedge from opposite side and protection of the papilla
- ▶ Provide effective isolation and reduce potential for trauma

## Wedge Refills

**M1008** Small (100)

**M1009** Medium (100)

**M1010** Large (100)



## WedgeShield (3 sizes)

- ▶ Ultra thin medical grade stainless steel panel
- ▶ Guard detaches after prep, leaving the wedge in place
- ▶ Facilitates faster and easier cutting of cavity and lowers risk of interproximal damage to adjacent tooth

## WedgeShield Refills

**M1011** Small (50)

**M1012** Medium (50)

**M1013** Large (50)




## Ring Forceps

- ▶ Innovative built-in width limiter to minimize ring damage
- ▶ The tips have a special angle with an occlusal-gingival orientation

## Instrument Refills

**M1015** Forceps Refill (1)



## Pin Tweezers

- ▶ Gold plated pin positively grips Smart-Loop™ in matrix and wedge for maximum control
- ▶ Default closed position prevent dropping matrix or wedge
- ▶ Convenient ball burnisher built-in to the tip

## Instrument Refills

**M1016** Pin Tweezers Refill (1)



7-145 Select Avenue,  
Toronto, ON, M1V 5M8  
Canada

Tel: 1-416-332-2222  
orders@marksondental.ca  
www.marksondental.ca



The specifications and appearance of the products are subject to change without notice. Printed colour can be slightly different from actual colour. Changes, variations and errors reserved.

(M1001FUR01)