

Your Dental Matrix Specialist





## **Flowdam**<sup>™</sup> Flowable Dental Dam

Flowdam is light-cured material barrier is a passively adhesive or sealing. Flowdam is specially formulated in green to ensure a safe, unmistakable protection every time. It offers ideal thyxotropy, polymerization time, and final hardness to create a perfect shape for restoration matrix system.





## Major Application

- Support anterior matrix system where interdental wedge would risk deforming matrix.
- Protects soft tissue with impervious seal such whitening and sandblasting procedure
- Recreate the occlusal anatomy with the stamp technique

## Key Feature

- Fast light cure.
- Stays cool when cured for patient comfort.
- Stays where placed and does not run. Remains flexible after curing. Tear resistant. Easy to remove.
- Stacks on itself smoothly and evenly.
- Moisture friendly chemistry works well in oral environment.
- Adheres to models for making bleaching trays.
- Blocks out undercuts.
- Repairs tears in rubber dams.



## Order information

M1040 M1033 Flowdam<sup>™</sup> Flowable Dental Dam, 4 Syringe + 25 Tips Flowdam<sup>™</sup> Syringe Tips Refill, 20 GA, (100)





7-145 Select Avenue. Canada

Tel: 1-416-332-2222 Toronto, ON, M1V 5M8 orders@marksondental.ca www.marksondental.ca

