





## **Finish with long lasting performance**

Achieving successful dental restorations involves excellent finishing and polishing. SurePoli finishing and polishing discs make it easier than ever to perform the job, ensuring that your restorations achieve a final outcome that matches natural dentition.



Superior Polishing Disc System

FLEXIBILITY: The bendability of the disc enables effortless access to difficult-to-reach areas, including the posterior region, for effective polishing.

COLOR CODED: Easy, colour-coded system allows for quick identification of grit sequence. Discs are also available in two different diameter sizes depending on your application

SIMPLICITY: Is made of polyester coated with aluminum oxide particles. Its abrasive layer has no extra adhesive coating. Once you attach the disc to the mandrel, you are instantly prepared for use!



Disc Type	Indicated for
Coarse	GROSS REDUCTION AND CONTOURING  • To adjust primary geometric form.
Medium	FINISHING  • To remove composite excess, shape occlusal anatomy, lingual fissures, and secondary anatomy.
Fine	POLISHING To erase surface scratches, and reduce surface roughness
Super-Fine	HIGH GLOSS POLISHING  • To reduce surface roughness to high gloss

## Great for composite, glass ionomers, amalgams, semiprecious and precious metals

## Order information

M1052	Intro Kit, (30) each size, (1) Mandrel
M1053	Disc Refills, Course, 12.7mm, 1/2" (85)
M1054	Disc Refills, Medium, 12.7mm, 1/2" (85)
M1055	Disc Refills, Fine, 12.7mm, 1/2" (85)
M1056	Disc Refills, Super Fine, 12.7mm, 1/2" (85)
M1057	Disc Refills, Course, 9.5mm, 3/8" (85)
M1058	Disc Refills, Medium, 9.5mm, 3/8" (85)
M1059	Disc Refills, Fine, 9.5mm, 3/8" (85)
M1060	Disc Refills, Super Fine, 9.5mm, 3/8" (85)
M1061	Mandrel RA Refill (3)





145 Select Avenue, Toronto, ON, M1V 5M8 Canada

Tel: 1-416-332-2222 info@marksondental.ca www.marksondental.ca

